



GETTING STARTED



WITH YOUR REZZIMAX TUNER PRO II

rezzimax.com/getting-started

TUNER TIP



INTENSITY (+)

ALGORITHMS

INTENSITY (-)

POWER

SCRAPING WING

MICRO USB PORT



TONGS

Stimulate multiple points simultaneously

REZZIMAX PILLOW



ACCESSORY STRAP & BELT PACK



FOAM BRAIN TIP

For extra-sensitive areas



CLEAR TIP

Softer option for sensitive areas

TUNER TIP - *Pressure points or to attach Tuner accessories*

INTENSITY (+) - *Increase the level of intensity of the Tuner*

ALGORITHMS - *Press to cycle through four pre-set algorithms (Blue, Green, Red and White)*

INTENSITY (-) - *Decrease the level of intensity of the Tuner (also decreases algorithms by 50%)*

POWER - *Press to turn off/on*

SCRAPING WING - *Perform Instrument Assisted Soft Tissue Mobilization (IASTM)*

MICRO USB PORT - *Plug in the provided Micro-USB cable and AC adapter to charge your Tuner*



ALGORITHMS

The Tuner will give two distinct vibrations at the completion of a cycle

BLUE - Calm

Calms the nervous system. It replicates the waves of the ocean with a rise and fall pattern. (4.5 minutes)

Use for overworked muscles, everyday aches & pains and dissipating mental stress.

GREEN - Accelerated Calm

Offers the same calming pattern as Blue Mode, only the cycle is accelerated at 3X faster speed (1.5 minutes)

Use for overworked muscles, everyday aches & pains, and dissipating mental stress.

RED - Wake Up

The pattern is reversed from the calm modes and stimulates the nervous system. (3.5 minutes)

Use for extra energy, neuropathy, autoimmune and chronic pain.

WHITE - Relax

A meandering program that ebbs & flows with no repeating pattern. (∞)

Use for meditation and inducing a state of restfulness.

Visit rezzimax.com/rezzipes to watch our easy to learn tutorials for the most common issues.



Head & Neck



Elbow & Hand



Foot & Ankle



& More!



SCAN NOW



REZZIPET FOR ANIMALS

rezzipet.com



Even more Rezzipes and resources made specifically for our four-legged friends!